ENGLISH

INSTRUCTIONS FOR USE AKILINK MINI INVASIVE ACHILLES TENDON REPAIR SYSTEM STERILE INSTRUMENTATION FOR FOOT SURGERY . SINGLE USE

In accordance with the directive 93/42/EEC relative to medical devices and its amendments, this product must be handled by WELL-TRAINED and QUALIFIED PERSONS, AWARE of these INSTRUCTIONS FOR USE.

1. Description of the medical device

The AKILINK mini invasive Achilles tendon repair system is a method to treat acute Achilles tendon ruptures. It is a minimally invasive procedure that allows direct visual control of the repair, as well as percutaneous introduction of the sutures. The surrounding soft tissues and tendon itself are always treated with the utmost care to avoid any local trauma. The AKILINK mini invasive Achilles tendon repair system includes:

- One AKILINK instrument,
- One needle pusher,
- One lateral extension,
- Two needles.

These medical devices are made out of polycarbonate, IXEF GS-1022 and stainless steel 316L in accordance with ASTM F138-139 and ISO5832-1 Standards. These medical devices are sold sterile. Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation. These devices do not contain phthalates unless this is indicated on the label.

2. Indications

The AKILINK system is indicated for the mini invasive repair of acute (less than 10 days) closed ruptures of the Achilles tendon, located between 2 cm and 8 cm above the tuberosity of the calcaneum.

3. Contraindications

The device should not be used in a patient who has currently, or who has history of:

- Acute or chronic, systemic inflammations,
- Active infections,
- Chronic rupture,
- Open rupture,
- Complex open rupture with skin defect,
- Rupture located between 0 and 2 cm above the tuberosity of the calcaneum or higher than 8 cm above the tuberosity of the calcaneum,
- Previous local surgery,
- Sensitivity/allergies to the medical device materials (cf paragraph 1),
- Patient under steroid, quinolones,
- Diabetic patient.

4. Warnings and precautions

Physician must determine if the use of the instrumentation is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition,
- Use of chemotherapy, etc.,
- Drug and/or alcohol and/or smoke addiction and/or abuse,
- Obesity,
- Compromised wound healing,
- Vascular disorder,
- A patient unwilling or unable to comply with postoperative instructions,
- Paediatric age.

5. Complications

Complications may include but are not limited to:

- Adhesions,
- Risk of re-rupture in particular in case of non-compliance to

postoperative instruction or postoperative trauma,

- Sural nerve injury,
- Skin necrosis,
- Pain, discomfort or abnormal sensations,
- Infections, hematoma, allergy, thrombosis, pulmonary embolism.

6. Use of the instrumentation

Knowledge of surgical techniques, proper reduction, and post-operative patient management are considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon. Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

6.1 Preoperative

- Read carefully the Surgical Technique.
- Check packaging and labelling integrity before use.
- Only the instruments designed and manufactured by In2Bones should be used.
- The sterility is guaranteed as long as the packing has not been damaged and before the end of the sterility validity.
- Do not use any sterile device for which the packaging has been opened or damaged outside the operating theatre. Inner packaging should be handled under sterile conditions (persons/instruments)

INFORMATION
Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer.

Last revision: 10/2022

Please note that a single use device which comes into contact with human blood or tissue should not be re-used and should be returned to the manufacturer or properly disposed. The company declines all responsibility in the event of such re-use.

Re-sterilization of devices sold sterile is forbidden.

7. Storage

Store in dry place.

8. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot under any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

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