ENGLISH INSTRUCTIONS FOR USE • DUAFIT® STERILE IMPLANTS FOR FOOT SURGERY - SINGLE USE

This product must be handled and/ or implanted by WELL-TRAINED and QUALIFIED PERSONS, AWARE of these INSTRUCTIONS FOR USE.

1. Description of the medical device

The DUAFIT® implant for arthrodesis is designed for proximal interphalangeal joint and is available in different sizes.

They are made out of PEEK (Poly Ether Ether Ketone) according to standard ASTM F2026. These medical devices are sold sterile.

Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation. These devices do not contain phtalates unless this is indicated on the label.

2. Indications

The DUAFIT® implant is intended for fixation of proximal interphalangeal joint arthrodesis of the lesser toes.

- Examples include:
- rigid or semi-rigid hammetoe deformity

- revision of failed arthroplasty or arthrodesis
- 2nd toe shortening.

3. Contraindications

The implant should not be used in a patient who has currently, or who has history of:

- acute or chronic inflammations, whether local or systemic,
- active infections,
- sensitivity/allergies to the implant materials (cf paragraph 1).

4. Warnings and precautions

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition
- Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic developments, acute osteopenia, bone tumor, etc...)
- Use of steroid derivatives, chemotherapy, ... etc.
- Drug and/or alcohol and/or smoke addiction and/or abuse
- Obesity
- Compromised wound healing
- Vascular disorder
- A patient unwilling or unable to comply with postoperative instructions

5. Complications

Complications may include but are not limited to:

- Pain, discomfort or abnormal sensations due to presence of the implant
- Bending, loosening and/or breakage
- Migration of the implant position
- Bone loss due to stress shielding
- Deformation recurrence, loss of correction
- Delayed union or pseudarthrosis
- Tissue irritation (tendon, nerve, soft tissue)
- Infections, hematoma, allergy, thrombosis

6. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon. Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

6.1. Preoperative

- Read carefully the Surgical Technique.
- Use x-rays to ensure that implant is able to be adapted to the joint concerned (corticals too thin, a diaphysis very wide or too narrow).

- Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.

7. Storage
Store in dry place.

8. MRI/SCANNER

The DUAFIT® implant is manufactured in PEEK, which is a non-ferromagnetic material. The patient should be asked to systematically mention that he/she was implanted with a surgical implant in the foot.

The DUAFIT® implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the DUAFIT® implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

9. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot under any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician. INFORMATION Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer. Last revision: 10/2022

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician. INFORMATION Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer. Last revision: 10/2022

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