ENGLISH INSTRUCTIONS FOR USE • NEOFIEW® STERILE IMPLANTS FOR HAND SURGERY - SINGLE USE

This product must be handled and/ or implanted by WELL-TRAINED and QUALIFIED PERSONS, AWARE of these INSTRUCTIONS FOR USE.

1. Description of the medical device

System of plates and screws for distal sizes. Implants existing in different sizes.

The NEOFIEW® plates have to be fixed with the NEO screws. The NEO screws exist in Locking (L) and Non-Locking (NL) versions.

The NEOFIEW plate shall be fixed with the NEO Locking (L) screws preferentially, and with a maximum of two (2) NEO Non-Locking screws on the whole construct.

The NEOFIEW® plates are made out of PEEK (Poly Ether Ether Ketone) according to standard ASTM F2026. The NEOFIEW® plates include radiopaque markers made out of Tantalum according to ASTM F600 for radiologic al evaluation. The NEOFIEW® plate is sold preassembled with a single use instrument (distal block) in its packaging.

The NEO screws are made out of titanium alloy Ti-6Al-4V according to standards ISO 5832-3 and ASTM F136.

The medical devices are sold sterile. Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation. These devices do not contain phthalates unless this is indicated on the label.

2. Indications

The NEOFIEW® Plate is intended for fixation of intra-articular and extra-articular fractures of the distal radius.

3. Contraindications

The implant should not be used in a patient who has currently, or who has history of:

- acute or chronic inflammations, whether local or systemic,
- active infections,
- sensitivity/allergies to the implant materials (cf paragraph 1).

4. Warnings and precautions

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition
- Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic degeneration, acute osteopenia, bone tumor, etc…)
- Use of steroid derivatives, chemotherapy,… etc.
- Drug and/or alcohol and/or smoke addiction and/or abuse
- Obesity
- Compromised wound healing
- Vascular disorder
- A patient unwilling or unable to comply with postoperative instructions

Complications may include but are not limited to:

- Pain, discomfort or abnormal sensations due to presence of the implant
- Bending, loosening and/or breakage
- Migration of the implant position
- Bone loss due to stress shielding
- Deformation recurrence, loss of correction
- Delayed union or pseudarthrosis
- Infections, hematoma, allergy, thrombosis,
- Tendon rupture, including late rupture,
- Tendon irritation/injuries,
- Heterotopic ossification,
- Compartment syndrome,
- Nerve irritation/injuries.

5. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon. Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

5.1. Preoperative

Read carefully the Surgical Technique.

Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.

Prepare all implants and instruments immediately for the surgical procedure. Do not attempt a surgical

procedure with non-functional, broken, suspect or damage instrument.

-The sterility is guaranteed as long as the packing has not been damaged and before the end of the sterility validity.

-The size and number of implants should be assessed based on the preoperative X-rays.

-After measurements, some implants from each adequate size should be made available to have a sufficient inventory for the surgery.

-An additional implant from each size should be made available to replace any implant that might be accidentally contaminated during the surgery.

-Before the first surgery, the surgeon and assistants should manipulate the instruments to familiarize themselves with the material.

5.2. Peroperative

-The surgery should be performed by a surgeon with adequate background in orthopaedics and with respect to the different steps described in the Surgical Technique.

-Implants should be handled with care to avoid any scratch (risk of incipient break).

-Under no circumstances should the implant be modified.

-Alternate fixation methods should be available intraoperatively.

-Opening of the instruments set must be done according to aseptic condition.

-Check packaging and labeling integrity before use.

-Do not use any implant for which the packaging has been opened or damaged outside the operating theatre. Inner packaging should be handled under sterile conditions (persons/instruments).

-Use an X-ray reader screen to verify the implants positioning.

-In order to prevent tendon rupture, the NEOFIEW® Plate should be positioned proximal to the watershed line.

-The NEOFIEW plate shall be fixed with the NEO Locking (L) screws preferentially, and with a maximum of two (2) NEO Non-Locking screws on the whole construct.

-Acute or chronic inflammations, whether local or systemic,

- active infections,

- sensitivity/allergies to the implant materials (cf paragraph 1).

4. Warnings and precautions

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition
- Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic degeneration, acute osteopenia, bone tumor, etc…)
- Use of steroid derivatives, chemotherapy,… etc.
- Drug and/or alcohol and/or smoke addiction and/or abuse
- Obesity
- Compromised wound healing
- Vascular disorder
- A patient unwilling or unable to comply with postoperative instructions

Complications may include but are not limited to:

- Pain, discomfort or abnormal sensations due to presence of the implant
- Bending, loosening and/or breakage
- Migration of the implant position
- Bone loss due to stress shielding
- Deformation recurrence, loss of correction
- Delayed union or pseudarthrosis
- Infections, hematoma, allergy, thrombosis,
- Tendon rupture, including late rupture,
- Tendon irritation/injuries,
- Heterotopic ossification,
- Compartment syndrome,
- Nerve irritation/injuries.

5. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon. Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

5.1. Preoperative

Read carefully the Surgical Technique.

Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.

Prepare all implants and instruments immediately for the surgical procedure. Do not attempt a surgical

removal should be followed by careful postoperative management to avoid re-fractures. If the patient is older and has a low activity level, the surgeon may elect not to remove the implant in order to eliminate the risks of another surgery.

7. Storage

Store in dry place.

8. MRI/SCANNER

The patient should be asked to systematically mention that he/she was implanted with a surgical implant.

The NEOFIEW implant has not been evaluated for safety and compatibility with the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the NEOFIEW implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

9. Information on the products / Responsibility:

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot under any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products.

The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

INFORMATION
Should any information regarding the products or their uses be required, please contact the representative distributor or distributor or directly contact the manufacturer.
Last revision: 10/2022

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