**ENGLISH**

**INSTRUCTIONS FOR USE • NeoFit® MTP plate system STERILE IMPLANTS FOR FOOT SURGERY - SINGLE USE**

This product must be handled and/ or implanted by WELL-TRAINED and QUALIFIED PERSONS, AWARE of these INSTRUCTIONS FOR USE.

## 1. Description of the medical device

Screw – plate system existing in different models, diameters and lengths. These medical devices are made out of Titanium Alloy Ti-6Al-4V in accordance with ASTM F136 and ISO5832-3 Standards.

The NeoFit® MTP plate has to be fixed with NeoFit® MTP screws available in locking or non-locking version. These medical devices are sold sterile. Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation. These devices do not contain phthalates unless this is indicated on the label.

### 2. Indications

The NeoFit® MTP plate system is indicated in the fixation of arthrosis of the distal metatarsophalangeal articulation, and in particular in the case of:

- Hallux rigidus
- Severe hallux valgus (IM angle >20° HV angle >40°)
- Deformation due to rheumatoid arthritis
- Post-trauma arthrosis
- Neuromuscular instability
- The addition of an oblique compression screw through the joint is necessary (for example IBS™ screw) in order to ensure the closure of the plantar hinge and reinforce the stability of the arthrodese.

### 3. Contraindications

The implant should not be used in a patient who has currently, or who has history of:

- acute or chronic, systemic inflammations,
- active infections,
- sensitivity/allergies to the implant materials (cf paragraph 7).

## 4. Warnings and precautions

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition
- Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic degeneration, acute osteopenia, bone tumor,...)
- Use of steroid derivatives, chemotherapy...
- Drug and/or alcohol and/or smoke addition and/or abuse
- Obesity
- Compromised wound healing
- Vascular disorder
- A patient unwilling or unable to comply with postoperative instructions
- Complications may include but are not limited to:
  - Pain, discomfort or abnormal sensations due to presence of the implant
  - Bending, loosening and/or breakage
  - Migration of the implant position
  - Bone loss due to stress shielding
  - Deformation recurrence, loss of correction
  - Delayed union or pseudarthrosis
  - Infections, hematoma, allergy, thrombosis

### 5. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations for a successful outcome. Criteria for patient selection are the responsibility of the surgeon.

## 6. Removal of the implant after healing

Particularly in young active patients, implants may loosen, fracture, migrate, increase the risk of infection, cause pain, or stress shield bone – even after normal healing. The surgeon should consider the risks and benefits when deciding whether or not to remove an implant. Implant removal should be followed by careful postoperative management to avoid re-fractures. If the patient is older and has a low activity level, the surgeon may elect not to remove the implant in order to eliminate the risks of another surgery.

Removal of In2Bones implant can be performed using standard ablation instrumentation and for screws removal, using torx screwdrivers provided for implantation.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

## 7. Storage

Store in dry place.

## 8. MRI Safety Information

Non-clinical testing and magnet resonance imaging (MRI) conditions were performed to evaluate the orthopedic implants from In2Bones.

contaminated during the surgery. Before the first surgery, the surgeon and assistants should manipulate the instruments to familiarize themselves with the material.

### 5.2. Preoperative

The surgery should be performed by a surgeon with adequate background in orthopaedics and with respect to the different steps described in the Surgical Technique.

Implants should be handled with care to avoid any scratch (risk of incipient break).

Under no circumstances should the implant be modified.

Alternate fixation methods should be available intraoperatively. Opening of the instruments set must be done according to aseptic condition.

Check packaging and labeling integrity before use.

Do not use any implant for which the packaging has been opened or damaged outside the operating theatre. Inner packaging should be handled under sterile conditions (persons/instruments).

Use laser marked indications on the plate to adequately position the system.

Do not apply excessive plate bending (no bending further than the indicated on associated bending instrument). Do not apply reverse, repeated and/or multidirectional bending.

Only bending instruments provided by In2Bones should be used to apply bending to In2Bones plates.

### 5.3. Postoperative

A rigid sole postoperative shoe should be used during 6 weeks.

The patient should be advised that a second more minor procedure for the removal of the implants may be necessary.

Some X-rays should be periodically done to check the postoperative progress and prevent any complication.

Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of excessive or early weight-bearing or excessive muscular activity. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until stability is evident.

The patient should be encouraged to report to his/her surgeon any unusual changes of the operated extremity. If evidence suggests loosening of the implant (particular pain and progressive changes in the radiographs) an intensified schedule of check-ups is advised and new warning and instructions to the patient may be appropriate regarding further activity restrictions.

The patient should be encouraged to receive prompt medical attention for any infection that could occur, such as the following: member level or elsewhere in the body.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Furthermore, re-use of single use device may create risks of contamination from one patient to another or the user. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even if the device appears undamaged, the implants may have defects and internal stress patterns that may cause material fatigue.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Furthermore, re-use of single use device may create risks of contamination from one patient to another or the user. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even if the device appears undamaged, the implants may have defects and internal stress patterns that may cause material fatigue.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Furthermore, re-use of single use device may create risks of contamination from one patient to another or the user. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even if the device appears undamaged, the implants may have defects and internal stress patterns that may cause material fatigue.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

The company declines all responsibility in the event of such re-use. Re-sterilization of devices sold sterile is forbidden.

## 9. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot undertake any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

WARRANTING (USA) This device has not been approved for fixation or maintaining purposes in the back part of the intervertebral body (neck of vertebra) or at cervical, thoracic or lumbar levels. CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

Non-clinical testing and electromag-netic simulations demonstrated that the orthopedic implants from In2Bones SAS are MR Conditional. A patient with one of these implants can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only

- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)

- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, an orthopedic implant from In2Bones is expected to produce a maximum temperature rise of less than or equal to 7°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the presence of implants from In2Bones extends approximately 10-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system

## 9. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot undertake any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

WARRANTING (USA) This device has not been approved for fixation or maintaining purposes in the back part of the intervertebral body (neck of vertebra) or at cervical, thoracic or lumbar levels. CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

### INFORMATION

Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer.

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