## ENGLISH

### INSTRUCTIONS FOR USE

#### STERILE KLIP™ KWIRE TIP PROTECTOR SINGLE USE

**This product must be handled by WELL-TRAINED and QUALIFIED PERSONS, AWARE of these INSTRUCTIONS FOR USE.**

#### 1. Description of the medical device

The medical devices covered by these INSTRUCTIONS FOR USE are sterile single use KLIP™ Kwire tip protectors.

KLIP™ Kwire tip protectors are external devices designed to protect the ends of the fingers or toes when realignment wire is used by the surgeon.

The Kwire tip protectors are available in two different sizes characterized by a color code:

- KLIP™ Kwire tip protector for wires diameter 0.8mm to 1.0mm : Orange color, identified by the letter “O”

- KLIP™ Kwire tip protector for wires diameter 1.2mm to 1.6mm : Violet color, identified by the letter “V”

The KLIP™ Kwire tip protectors are made out of high density polyethylene. They may be supplied either individually packaged, or packed by several units (example pack of 5 units, or 10 units).

These medical devices are sold sterile. Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation.

#### 2. Indications

KLIP™ Kwire tip protector is indicated after insertion of realignment wire with diameter ranging from 0.8mm to 1.6mm. KLIP™ Kwire tip protector is intended to protect per-operatively and post-operatively the patient and the staff from sharp ends and snagging wire. KLIP™ Kwire tip protector also prevents a deeper insertion of the wire into the patient.

KLIP™ Kwire tip protector facilitates the removal of the wire.

These medical devices should not be implanted.

#### 3. Use of the medical device

Only medical professionals who are thoroughly familiar with the function, application, and use should use this medical device. Check packaging and labeling integrity before use. The sterility is guaranteed as long as the packaging has not been damaged (film scratched, label missing, questionable packing...) and before the expiration date. Do not use any product for which the packaging has been opened or damaged. KLIP™ Kwire tip protector is intended for use per or post-operatively following wire insertion into the patient. The size of the KLIP™ Kwire tip protector must be chosen according to the diameter of the wire implanted into the patient:

- For wire diameter from 0.8mm to 1.0mm: Use the orange KLIP™ Kwire tip protector.

- For wire diameter from 1.2 mm to 1.6mm: Use the violet KLIP™ Kwire tip

protector.

The wire must be bent with an angle of 90°. Leave a minimum distance of 12mm between the patient and the place of bending. The KLIP™ Kwire tip protector is inserted and locked on the wire according to the following surgical technique:

**1 : Slide the KLIP™ Kwire tip protector on the wire under the bridge**

**2 : Rotate the KLIP™ Kwire tip protector around the bent axis of the wire**

**3 : Fold the two sides of the KLIP™ Kwire tip protector on each another until clipping**

For the removal of the wire, execute first a rotation of the Kwire tip protector in order to release the wire from bone and soft tissues. Then pull the Kwire tip protector in order to remove the wire. Once the wire is removed, the KLIP™ Kwire tip protector must be disposed of in accordance to hospital protocol.

#### 4. Re-sterilization / Re-use

Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure that may result in patient injury or illness. Furthermore, re-use of single use devices may create risks of contamination from one patient to another or the user. Please note that a single use device which comes into contact with human blood or tissue should not be re-used and should be returned to the manufacturer or properly disposed. Even though they may appear undamaged, such devices may have defects and internal stress patterns that may cause material fatigue. The company declines all responsibility in the event of such re-use.

#### 5. Storage

Store in dry place.

#### 6. Product information disclosure / Liability

In2Bones® has exercised reasonable care in the selection of materials and the manufacture of these products and warrants that

the products are free from manufacturing defects. In- aBones® excludes all other warranties, whether expressed or implied, including but not limited to, any implied warranties of merchantability or fitness for a particular purpose. In2Bones® shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from use of this product.

In2Bones® neither assumes nor authorizes any person to assume for it any other or additional liability or responsibility in connection with these products. In2Bones® intends that this device should be used only by medical professionals who are thoroughly familiar with the function, application, and use of the device.

#### 7. MRI Safety Information to the attention of surgeons OUTSIDE of the USA only

**MR** The KLIP™ device from In2Bones® is MR Safe, which is an item that poses no known hazards in all MR environments.

**INFORMATION**

Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer. Last revision : 10/2022

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