ENGLISH INSTRUCTIONS FOR USE **NEOFUSE ANKLE FUSION PLATING SYSTEM** STERILE IMPLANTS FOR FOOT SURGERY - SINGLE USE

In accordance with the directive 93/42/EEC relative to medical devices and its amendments, this product must be handled and/or implanted by WELL-TRAINED and QUALIFIED PERSONS, AWARE OF THESE INSTRUCTIONS FOR USE.

1. Description of the medical device

System of plates and screws for ankle fusion.

The NEOFUSE Ankle Fusion plates have to be fixed with the associated :

- Talar fixation using fixed angle locking system (screw and cap) of 3.5mm diameter.

- Tibial fixation using fixed angle locking system (screw and cap) of 4.5mm diameter

- 4.5mm diameter transfoal screw.

These medical devices are made out of Titanium Alloy Ti-6Al-4V in accordance with ASTM F136 and ISO5832-3 Standards.

These medical devices are ISO sterile. Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation.

These devices do not contain phtalates unless this is indicated on the label.

2. Indications

The NEOFUSE Ankle Fusion plating system is indicated for anterior fixation of ankle arthroses and fractures, including the distal tibia, talus and calcaneus.

The addition of a compression screw through the tibiotalar joint (for example IBSTM 6.5mm screw) is necessary in order to ensure the closure of the posterior hinge and reinforce the stability of the arthrodesis.

3. Contraindications

The implant should not be used in a patient who has currently, or who has history of:

- Acute or chronic, systemic inflammations,
- Active infections,
- Sensitivity/allergies to the implant materials (cf paragraph 1),
- Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic development, acute osteopenia, bone tumor, etc..)

4. Warnings and precautions

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Lacks good general physical condition
- Use of steroid derivatives, chemotherapies, etc.
- Drug and/or alcohol and/or smoke addiction and/or abuse
- Obesity
- Compromised wound healing
- Vascular disorder
- A patient unwilling or unable to comply with postoperative instructions

5. Complications

Complications may include but are not limited to:

- Pain, discomfort or abnormal sensations due to presence of the implant
- Bending, loosening and/or breakage
- Migration of the implant position
- Bone loss due to stress shielding
- Deformation recurrence, loss of correction
- Delayed union or pseudarthrosis
- Tissue irritation (tendon, nerve, soft tissue)
- Adjacent joint osteoarthritis
- Infections, hematoma, allergy, thrombosis

6. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon. Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

6.1. Preoperative

- Read carefully the Surgical Technique.
- Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.

- Prepare all implants and instruments necessary for the surgical procedure. Do not attempt a surgical procedure with non-functional, broken, suspect or damage instrument.

- The sterility is guaranteed as long as the packing has not been damaged and before the end of the sterility validity.

- The size and number of implants should be assessed based on the preoperative X-rays.

- After measurements, some implants from each adequate size should be made available to have a sufficient inventory for the surgery.

- An additional implant from each size should be made available to replace any implant that might be accidentally contaminated during the surgery.

- Before the first surgery, the surgeon and assistants should manipulate the instruments to familiarize themselves with the material.

- In case of severe deformity or bone defects, reconstruction with allograft or autograft from iliac crest or other anatomic regions may be necessary.

- Excessive length of the fibula causing lateral impingement may make shortening necessary. This can be easily done through the same anterior approach.

6.2. Preoperative

- The surgery should be performed by a surgeon with adequate background in orthopaedics and with respect to the different steps described in the Surgical Technique.

- Implants should be handled with care to avoid any scratch (risk of incision break).

- Under no circumstances should the implant be modified.

- Alternate fixation methods should be available intraoperatively.

- Opening of the instruments set must be done according to aseptic condition.

- Check packaging and labeling integrity before use.

- Do not use any implant for which the packaging has been opened or damaged outside the operating theatre. Inner packaging should be handled under sterile conditions (persons/instruments).

- Do not apply any plate bending.

- The NEOFUSE Ankle Fusion plates have to be fixed with the associated :

- Talar fixation using fixed angle locking system (screw and cap) of 3.5mm diameter.
- Tibial fixation using fixed angle locking system (screw and cap) of 4.5mm diameter
- 4.5mm diameter transfoal screw.

No other fixation screw should be used.

6.3. Postoperative

- The patient should be advised that a second more minor procedure for the removal of the implants may be necessary.

- Some X-rays should be periodically done to check the postoperative progress and prevent any complication.

- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of excessive or early weight-bearing or excessive muscular activity. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until stability is evident.

- The patient should be encouraged to report to his/her surgeon any unusual changes of the operated extremity. If evidence suggests loosening of the implant (particular pain and progressive changes in the radiographs) an intensified schedule of check-ups is advised and new warning and instructions to the patient may be appropriate regarding further activity restrictions.

- The patient should be encouraged to receive prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.
- The patient should be encouraged to report to prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.
- Adjacent joint osteoarthritis
- Infections, hematoma, allergy, thrombosis

6.4. Re-use / Re-sterilization
Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Furthermore, re-use of single use devices may create risks of contamination from one patient to another or the user. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even though they may appear undamaged, the implants may have been contaminated by internal stress patterns that may cause material fatigue. The company declines all responsibility in the event of such re-use.

Re-sterilization of devices sold sterile is forbidden.

- Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.

- Prepare all implants and instruments necessary for the surgical procedure. Do not attempt a surgical procedure with non-functional, broken, suspect or damage instrument.

- The sterility is guaranteed as long as the packing has not been damaged and before the end of the sterility validity.

- The size and number of implants should be assessed based on the preoperative X-rays.

- After measurements, some implants from each adequate size should be made available to have a sufficient inventory for the surgery.

- An additional implant from each size should be made available to replace any implant that might be accidentally contaminated during the surgery.

- Before the first surgery, the surgeon and assistants should manipulate the instruments to familiarize themselves with the material.

implant. Implant removal should be followed by careful postoperative management to avoid re-fractures. If the patient is older and has a low activity level, the surgeon may elect not to remove the implant in order to eliminate the risks of another surgery.

8. Storage

Store in dry place.

9. MRI/SCANNER

The patient should be asked to systematically mention that he/she was implanted with a surgical implant in the foot.

The NEOFUSE Ankle Fusion Plating System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the NEOFUSE Ankle Fusion Plating System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

10. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal or clinical, thoracic or lumbar levels. implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot under any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products. The intention of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

WARNING (USA): This device has not been approved for fixation or maintaining purposes in the back part of the intervertebral body (neck or vertebra) or at cervical, thoracic or lumbar levels.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

- Talar fixation using fixed angle locking system (screw and cap) of 3.5mm diameter.

- Tibial fixation using fixed angle locking system (screw and cap) of 4.5mm diameter

- 4.5mm diameter transfoal screw.

No other fixation screw should be used.

6.3. Postoperative

- The patient should be advised that a second more minor procedure for the removal of the implants may be necessary.

- Some X-rays should be periodically done to check the postoperative progress and prevent any complication.

- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of excessive or early weight-bearing or excessive muscular activity. Postoperative care and physical therapy should be structured to prevent loading of the operative extremity until stability is evident.

- The patient should be encouraged to report to his/her surgeon any unusual changes of the operated extremity. If evidence suggests loosening of the implant (particular pain and progressive changes in the radiographs) an intensified schedule of check-ups is advised and new warning and instructions to the patient may be appropriate regarding further activity restrictions.

- The patient should be encouraged to receive prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.
- The patient should be encouraged to report to prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.
- Adjacent joint osteoarthritis
- Infections, hematoma, allergy, thrombosis

6.4. Re-use / Re-sterilization
Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death. Furthermore, re-use of single use devices may create risks of contamination from one patient to another or the user. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even though they may appear undamaged, the implants may have been contaminated by internal stress patterns that may cause material fatigue. The company declines all responsibility in the event of such re-use.

Re-sterilization of devices sold sterile is forbidden.

- Only the instruments designed and manufactured by In2Bones should be used in combination with the implants.

- Prepare all implants and instruments necessary for the surgical procedure. Do not attempt a surgical procedure with non-functional, broken, suspect or damage instrument.

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