ENGLISH

MRI Safety Information to the attention of surgeons OUTSIDE of the USA only

MR Conditional

Non-clinical testing and electromagnetic simulations were performed to evaluate the orthopedic implants from In2Bones® (Table 1). Non-clinical testing and electromagnetic simulations demonstrated that the orthopedic implants from In2Bones® are MR Conditional. A patient with one of these implants can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, an orthopedic implant from In2Bones® is expected to produce a maximum temperature rise of less than or equal to 7.0°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the orthopedic implants from In2Bones® extends approximately 10-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

Table 1.

NeoFit MTP Plate System
IBS® Screw System
TriWay TTC Fusion Nail
Anterior Ankle Fusion Plate System
NeoLink Plating System
NeoView Plating System

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