ENGLISH

INSTRUCTIONS FOR USE • QUANTUM® TAR STERILE IMPLANTS FOR ANKLE SURGERY • SINGLE USE

These products must be handled and/or implanted by WELL-TRAINED and QUALIFIED PERSONS. AWARE OF THESE INSTRUCTIONS FOR USE.

1. Description of the medical device

The QUANTUM® total ankle prosthesis is a fixed-bearing, semi-constrained ankle prothesis with 2 components and composed of:

- a titanium component of association of a metallic (Ti6AlV) tibial tray implant fixed to a polymer (UHMWPE) insert,
- a metallic (CoCr) tarl implant reproducing the talus dome anatomy.

When all components are implanted, with the insert rigidly locked to the tibial tray, the polyethylene insert acts as a single bearing along the talor dome, enabling flexion-extension and rotation movement at the replaced ankle joint.

Each component of the QUANTUM® total ankle prothesis exists in different sizes and models. The minimal thickness available for the QUANTUM® fixed bearing insert is 5mm.

Recommended association of sizes of each component of the QUANTUM® total ankle prosthesis are described in the figure below:

TIBIAL TRAY > SIZE X OR SIZE X+1

TIBIAL INSERT > SIZE X

TALAR IMPLANT > SIZE X

Caution: the associated QUANTUM® talar implant and insert shall be of the same size. The QUANTUM® Tibial tray shall be of same size or one size upper than the talar implant.

The QUANTUM® tibial tray is manufactured out of titanium alloy Ti6Al4V compliant with ISO 5832-3 & ASTM F136 standards with porous Titanium coating according to ISO 13179-1. The QUANTUM® fixed insert implant is manufactured out of Ultra High Molecular Weight Polyethylene (UHMWPE) compliant with ISO 5834-2 & ASTM F648 standard.

The QUANTUM® talar implant is manufactured out of titanium chrome alloy (CoCr) compliant with ISO 5832-4 & ASTM F75 standards with porous Ttanium coating according to ISO 13179-1. These medical devices are sold sterile.

Elements sterilized using irradiation have been exposed to a minimum of 25kGy of gamma irradiation.

These devices do not contain phthalates unless this is indicated on the label.

2. Indications

The QUANTUM® total ankle prosthesis is indicated for total ankle replacement in primary or revision surgery for patients with ankle joints damaged by severe rheumatoid, post-traumatic, or degenerative arthritis. Note: In the United States, the ankle prosthesis is intended for cement use only.

3. Contraindications

The QUANTUM® Total Ankle Prosthesis is contraindicated for the following conditions:

- Sepsis,
- Active /prior deep infection in ankle/joint or adjacent bones
- Fever and/or local inflammation,
- Avascular necrosis of the talus / tibia,
- Osteoporosis / osteopenia,
- Poor skin coverage / soft-tissue quality around the ankle joint that would make the procedure unjustifiable,
- Inadequate or insufficient quality of bone stock
- Impaired joint laxity or tendon dysfunction,
- Neuromuscular or mental disorders which might jeopardize fixation and post-operative care,
- The patient must be advised of the limitations of the reconstruction and the need for protection of the prosthesis from full weight bearing until adequate fixation and healing have occurred.
- Excessive activity and trauma affecting the joint replacement have been implicated with failure of the reconstruction by loosening, fracture and wear of the prosthetic components. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.
- The patient should be cautioned to limit activities and protect the replaced joint from unreasonable stresses, and follow the instructions of the physician with respect to follow-up care and treatment. The patient should be closely monitored if a change at the operative site has been detected. The possibility of deterioration of the joint should be evaluated and possible revision surgery considered.
- The patient should be warned of surgical risks, and made aware of possible adverse effects. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can not be heavily damaged as a result of certain activity or trauma, has a finite expected service life, and may need to be replaced at some time in the future. The patient should also be advised of other risks that the surgeon believes should be disclosed.
- The patient should be advised that any noise or unusual sensation should be reported to the surgeon as it may indicate implant malfunction.
- Periodic post-operative x-rays are recommended for close comparison with early post-posed conditions to detect long term evidence of changes in position, loosening, bending, or cracking of components.
- The patient should be encouraged to receive prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.

4. Warnings and precautions

The following conditions tend to adversely affect the performance of the device:

- Obesity or excessive patient weight,
- Manual labor,
- Active sports participation and/or high activity level,
- Severe malalignment or instability that is not surgically correctable
- Alcohol and/or drug addiction and / or abuse,
- Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g. diabetes, steroid usage, immunosuppressive treatments).
- Severe deformities of the joint,
- Tumors of the supporting bone structures,
- Sensitivity, allergy or other reactions to implant materials (i.e. polyethylene, bone cement or metal),

5. Complications

- Malalignment/Incorrec component positioning,
- Intraoperative or postoperative bone fracture,
- Tissue irritation, nerve, tendon, soft tissue),
- Nerve injury,
- Active sports participation and/or high activity level,
- Severe malalignment or instability that is not surgically correctable
- Alcohol and/or drug addiction and / or abuse,
- Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g. diabetes, steroid usage, immunosuppressive treatments).
- Severe deformities of the joint,
- Tumors of the supporting bone structures,
- Sensitivity, allergy or other reactions to implant materials (i.e. polyethylene, bone cement or metal),

- Secondary necrosis of the talus / tibia,
- Other diseases affecting the quality of the soft tissue,
- Sensitivity, Allergy or Other reactions to implant materials, (i.e. polyethylene, bone cement or metal),
- Adjacent joint osteoarthritis.

6. Use of the implant

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management and considerations essential to a successful outcome. Criteria for patient selection are the responsibility of the surgeon.

Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

- 6.1. Preoperative**
- Read carefully the Surgical Technique.
- Only the instruments designed and manufactured by In2Bones must be used in combination with the implants.
- Prepare all implants and instruments necessary for the surgical procedure. Do not attempt a surgical procedure with non-functional, broken, suspect or damage instrument.
- The sterility is guaranteed as long as the packing has not been damaged and before the end of the sterility validity.
- Before the surgery the surgeon must check with the templates delivered with the system, the sizes of the implants that may be implemented during the surgery.
- After measurements, some implants from each adequate size (and possible adjacent sizes) should be made available to have a sufficient inventory for the surgery.
- An additional implant from each size should be made available to replace any implant that might be accidentally contaminated during the surgery.

Before the first surgery, the surgeon and assistants should receive appropriate training from In2Bones representative regarding the QUANTUM® total ankle prosthesis system to ensure thorough understanding of the implantation technique and the instrumentation, and familiarize themselves with the material.

- The surgeon shall check that the quality of the procedure is satisfactory enough to support the implantation.
- The surgery should be performed by a surgeon with adequate background in orthopaedics and with respect to the different steps described in the Surgical Technique.
- The correct selection of the type and size of the implant appropriate to the patient and the positioning of the implant are extremely important.
- Implants should be handled with care to avoid any scratch risk of incipient break or early wear). The prostheses must not be used if their functional surfaces have been damaged or have undergone shock, abrasion, or other deterioration. In case of revision, particular precautions should be taken not to damage the components that are not removed.
- Under no circumstances should the implant be modified.

The surgeon shall check that the quality of the bone is satisfactory enough to support the implantation.

- Alternate fixation methods should be available intraoperatively.
- Opening of the instruments set must be done according to aseptic conditions.
- Check packaging and labelling integrity before use.
- Do not use any implant for which the packaging has been opened or damaged outside the operating theatre. Inner packaging should be handled under sterile conditions (persons/instruments).

- 6.3. Postoperative**
- Detail instructions on the use and limitations of the device should be given to the patient.
- The patient must be advised of the limitations of the reconstruction and the need for protection of the prosthesis from full weight bearing until adequate fixation and healing have occurred.
- Excessive activity and trauma affecting the joint replacement have been implicated with failure of the reconstruction by loosening, fracture and wear of the prosthetic components. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.
- The patient should be cautioned to limit activities and protect the replaced joint from unreasonable stresses, and follow the instructions of the physician with respect to follow-up care and treatment. The patient should be closely monitored if a change at the operative site has been detected. The possibility of deterioration of the joint should be evaluated and possible revision surgery considered.
- The patient should be warned of surgical risks, and made aware of possible adverse effects. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can not be heavily damaged as a result of certain activity or trauma, has a finite expected service life, and may need to be replaced at some time in the future. The patient should also be advised of other risks that the surgeon believes should be disclosed.
- The patient should be advised that any noise or unusual sensation should be reported to the surgeon as it may indicate implant malfunction.
- Periodic post-operative x-rays are recommended for close comparison with early post-posed conditions to detect long term evidence of changes in position, loosening, bending, or cracking of components.
- The patient should be encouraged to receive prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.

6.4. Re-use / Re-sterilization

Products intended for single use must not be re-used (see symbols). Re-use may compromise the structural integrity of the device and/or lead to a failure of the device which may result in injury, illness or death. Furthermore, re-use of single use device may create risks of contamination from one technique chirurgicale et de l'instrumentation, et

patient to another or of the user. Any implant that has been soiled by blood, tissue, and/or fluid/matter should never be re-used. It must be handled according to hospital protocol. Even though they may appear undamaged, the implants may have defects and internal stress patterns that may cause material failure.

The company declines all responsibility in the event of such re-use.

Re-sterilization of devices sold sterile is forbidden.

7. Storage

Store in dry place.

8. MRI

The patient should be asked to systematically mention that he/she was implanted with a surgical implant in the ankle.

MR Conditional Non-clinical testing and electromagnetic simulations were performed to evaluate the orthopedic implants from In2Bones®. Nonclinical testing and electromagnetic simulations demonstrated that the orthopedic implants, including the QUANTUM® total ankle prosthesis, from In2Bones® are MR Conditional for use with all MRI implants can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5Tesla and 3Tesla, only.
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2W/kg (whole body) of scanning in the Normal Operating Mode.
- Under the scan conditions defined, an orthopedic implant from In2Bones® is expected to produce a maximum temperature rise of less than or equal to 7.0 °C after 15-minutes of continuous scanning. In non-clinical testing, the image contrast caused by the orthopedic implants from In2Bones® extends approximately 10mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

9. Information on the products / Responsibility

In2Bones has taken reasonable precautions in the selection of materials and in the manufacture of these products. However, In2Bones excludes any legal guarantee, whether express or implicit, including but not limited to, any implicit guarantee of the marketable quality or suitability for a specific use. In2Bones cannot under any circumstances be held responsible for any loss, damage or related costs or incidents, directly or indirectly linked to the use of this product.

In2Bones does not assume, and does not authorize any third party to assume on its behalf, any other responsibilities relating to these products. The selection of In2Bones is that this device should be used only by doctors having received appropriate training in techniques of orthopaedic surgery for its use.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

INFORMATION Detailed information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer. Last revision: 10/2022

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